



# Characterization of Coal and Its By-products Using Borate Fusions and ICP-OES Analyses

By **Marie-Ève Provencher, M.Sc. Chemist**  
USA Sales Manager



# Samples from Coal Industry



© Claisse, 2015. All Rights Reserved / Tous Droits Réservés.

Coal is mostly used by power plants to generate electricity

Emissions from power plants:

- Gas (fly ash)
- Solid waste (ash)

Why is trace analysis important?

- It helps to know the content of the coal by-products (coal ash and coal fly ash)
- It allows to know the content of the waste from power plants (coal ash and coal fly ash) for an eco-friendly disposal

## Certified Reference Materials (CRMs)

- 3 CRMs from coal industry
- High silica and alumina content
- Sample preparation method
  - Acid digestion (not efficient)
  - Harsh acid (HF) (not entirely safe)
  - Fusion

|           | Coal        | Coal ash   | Coal fly ash   |
|-----------|-------------|------------|----------------|
|           | NCS FC28127 | VS-7177-95 | EOP 12-1-02    |
|           | China       | Russia     | Czech Republic |
| <b>Al</b> | 3.47        | 14.3       | 16.1           |
| <b>Si</b> | 5.61        | 27.4       | 22.9           |

## Norms

### ASTM D6349

- Sample preparation:
  - Borate fusion
  - Acid digestion (HF/HNO<sub>3</sub>/HCl)
- ICP-OES
- **Al, Ba, Ca, Fe, K, Mg, Mn, Na, P, S, Si, Sr and Ti**

### AS 1038.14.1

- Sample preparation:
  - borate fusion
- AA
- **Al, Ba, Ca, Fe, K, Mg, Mn, Na, P, S, Si, Sr, Ti and Zn**

## Pretreatment of Samples

- Ramp up to 500 °C for 1 hour
- Ramp up to 750 °C for 2 hours
- Hold at 750 °C for 2 hours



| Coal        | Coal ash   | Coal fly ash   |
|-------------|------------|----------------|
| NCS FC28127 | VS-7177-95 | EOP 12-1-02    |
| China       | Russia     | Czech Republic |
| LOI         | 74.3%      | 0.42%          |
|             |            | 0.29%          |

# Sample Preparation by Borate Fusion

6

XRF Analysis

ICP and AA Analysis




# Sample Preparation by Borate Fusion

In a Pt/Au (95%/5%) crucible:

- 0.100 g of ashed sample
- 0.150 g LiNO<sub>3</sub>
- 1.000 g LiM/1.5%LiBr

Keep in mind:

- Low sample/flux ratio (1/10) for an optimized sample dissolution
- The non-wetting agent allows a complete pouring
- Silica and alumina are very soluble in lithium metaborate



# Fusion Technique Instruments

8

## 1. Manual

- Meeker burner and rising stand
- Muffle furnace
- Home made

## 2. Automatic

- Gas fired
- Electric



# Claisse® Automatic Fluxers

**ELECTRIC**



LeNeo® Claisse® fluxer



Eagon 2® Claisse® fluxer



TheOx® Claisse® fluxer

**GAS**



M4™ Claisse® fluxer

← Propane & Natural Gas

# Automated Fusion

- LeNeo® fusion instrument
- One fusion position
- Programmable parameters (temperature, duration and agitation):
  - Oxidation step at 800 °C
  - Fast temperature ramp-up (up-to 1200 °C)
  - The fusion steps are completed in 7 minutes



| Steps     | Duration (min) | Rocking angle | Rocking speed | Temperature (°C) |
|-----------|----------------|---------------|---------------|------------------|
| Oxidation | 2              | 0             | 0             | 800              |
| Fusion    | 1              | 0             | 0             | 1000             |
| Fusion    | 1              | 0             | 0             | 1050             |
| Fusion    | 1              | 20            | 15            | 1050             |
| Fusion    | 1              | 30            | 20            | 1050             |
| Fusion    | 1              | 40            | 30            | 1050             |
| Pouring   | -              |               |               |                  |
| Stirring  | 6              |               |               |                  |

# Automated Dissolution

The dissolution is completed in 6 minutes:

- Complete pouring
- 90 mL of 10% HNO<sub>3</sub>
- Magnetic agitation



11

| Steps     | Duration (min) | Rocking angle | Rocking speed | Temperature (°C) |
|-----------|----------------|---------------|---------------|------------------|
| Oxidation | 2              | 0             | 0             | 800              |
| Fusion    | 1              | 0             | 0             | 1000             |
| Fusion    | 1              | 0             | 0             | 1050             |
| Fusion    | 1              | 20            | 15            | 1050             |
| Fusion    | 1              | 30            | 20            | 1050             |
| Fusion    | 1              | 40            | 30            | 1050             |
| Pouring   | -              |               |               |                  |
| Stirring  | 6              |               |               |                  |

## Manual Dilution

- Complete to 100 mL
- Add the internal standard
- Dilution 1 in 2
- Final volume at 200 mL

**0.5% of total dissolved salt**



# Analysis by ICP-OES

## Instrumentation

- Perkin Elmer® Optima® 7300 DV ICP-OES

## Analysis

- Matrix-matched calibration curve
- Internal standards added



|                         |  |
|-------------------------|--|
| <b>Nebulizer</b>        | Gem Tip Cross flow   |
| <b>Spray chamber</b>    | Scott  |
| <b>Injector</b>         | Alumina 2 mm i.d.  |
| <b>RF</b>               | 1500 W   |
| <b>Argon flow</b>       | Plasma: 16 L/min<br>Nebulizer: 0.8 L/min<br>Auxiliary: 0.4 L/min |
| <b>Sample flow rate</b> | 1.0 mL/min   |

# Internal Standards

Internal standard associated with the element:

- Similar ionization potential
- Similar ionization state

| Elements | Elements of interest |                      |                    | Associated internal standards |                      |                  | Both         |  |
|----------|----------------------|----------------------|--------------------|-------------------------------|----------------------|------------------|--------------|--|
|          | Wavelength           | Ionization potential | Internal standards | Wavelength                    | Ionization potential | Ionization state | Viewing mode |  |
| Al       | 237.313              | 6.0                  | Ga                 | 294.364                       | 6.0                  | I                | Axial        |  |
| Ca       | 422.673              | 6.1                  | Eu                 | 381.967                       | 5.7                  | II               | Radial       |  |
| Fe       | 238.863              | 7.9                  | Re                 | 197.248                       | 7.9                  | II               | Axial        |  |
| K        | 766.490              | 4.3                  | Rb                 | 780.023                       | 4.2                  | I                | Axial        |  |
| Mg       | 279.077              | 7.6                  | Re                 | 197.248                       | 7.9                  | II               | Axial        |  |
| Na       | 589.592              | 5.1                  | Ga                 | 294.364                       | 6.0                  | I                | Radial       |  |
| Si       | 252.851              | 8.2                  | Te                 | 214.281                       | 9.0                  | I                | Axial        |  |
| Ti       | 337.279              | 6.8                  | Eu                 | 381.967                       | 5.7                  | II               | Axial        |  |

# Calibration, Method Detection Limits (MDLs) and Linearity

| Elements | Wavelength (nm) | View   | MDL (mg/L) | Calibration range |       | Correlation coefficient (pts = 5) |
|----------|-----------------|--------|------------|-------------------|-------|-----------------------------------|
|          |                 |        |            | (Min)             | (Max) |                                   |
| Al       | 237.313         | Axial  | 0.07       | 13.8              | 275   | 1.0000                            |
| Ba       | 413.065         | Axial  | 0.006      | 0.25              | 5.00  | 1.0000                            |
| Ca       | 422.673         | Radial | 0.1        | 6.00              | 115   | 1.0000                            |
| Fe       | 238.863         | Axial  | 0.02       | 6.25              | 125   | 1.0000                            |
| K        | 766.490         | Axial  | 0.002      | 0.900             | 18.0  | 0.9997                            |
| Mg       | 279.077         | Axial  | 0.007      | 0.880             | 17.5  | 1.0000                            |
| Mn       | 293.305         | Axial  | 0.001      | 0.138             | 2.76  | 1.0000                            |
| Na       | 589.592         | Radial | 0.05       | 0.280             | 5.50  | 1.0000                            |
| P        | 177.434         | Axial  | 0.03       | 1.60              | 30.0  | 0.9995                            |
| S        | 181.975         | Axial  | 0.06       | 2.80              | 55.0  | 1.0000                            |
| Si       | 252.851         | Radial | 0.07       | 15.0              | 300   | 1.0000                            |
| Sr       | 460.733         | Axial  | 0.002      | 0.330             | 6.50  | 0.9999                            |
| Ti       | 337.279         | Axial  | 0.002      | 2.80              | 55.0  | 1.0000                            |
| Zn       | 213.857         | Axial  | 0.003      | 0.150             | 3.00  | 1.0000                            |



## Results

16

- The method has been tested with 3 different types of samples:
  - Fly ash
  - Coal ash
  - Coal
- Accuracy and precision were calculated based on the certified reference value of each CRM (Al, Ca, Fe, K, Mg, Na, Si and Ti)
- Recovery was calculated with samples spiked before fusion (Ba, K, Mn, P, S, Sr and Zn)



# Precision and Accuracy in Brown Coal Fly Ash <sup>17</sup>

## EOP 12-1-02

| Elements | Certified values (%) | Experimental values (%) n=10 | MDL (%) | Precision $t_{0,975;9}$ (%) | Accuracy (%) | RSD (%) |
|----------|----------------------|------------------------------|---------|-----------------------------|--------------|---------|
| Al       | 16.1                 | 16.3                         | 0.01    | 0.1                         | 98.7         | 0.5     |
| Ca       | 1.49                 | 1.44                         | 0.02    | 0.04                        | 96.5         | 3.9     |
| Fe       | 5.17                 | 5.21                         | 0.003   | 0.02                        | 99.2         | 0.5     |
| K        | 0.651                | 0.630                        | 0.0004  | 0.011                       | 96.7         | 2.4     |
| Mg       | 0.581                | 0.623                        | 0.001   | 0.002                       | 93.0         | 0.5     |
| Na       | 0.361                | 0.37                         | 0.01    | 0.01                        | 97.6         | 3.5     |
| Si       | 22.9                 | 23.8                         | 0.01    | 0.4                         | 96.3         | 3.7     |
| Ti       | 3.61                 | 3.44                         | 0.0003  | 0.01                        | 95.3         | 0.5     |

Dilution factor: 1994  
n = 10

# Precision and Accuracy in Coal

18

| NCS FC28127 |                      |                              |         |                             |              |         |
|-------------|----------------------|------------------------------|---------|-----------------------------|--------------|---------|
| Elements    | Certified values (%) | Experimental values (%) n=10 | MDL (%) | Precision $t_{0,975;9}$ (%) | Accuracy (%) | RSD (%) |
| Al          | 3.47                 | 3.37                         | 0.003   | 0.03                        | 97.1         | 1.2     |
| Ca          | 1.88                 | 1.89                         | 0.006   | 0.04                        | 99.3         | 2.8     |
| Fe          | 1.02                 | 1.05                         | 0.0008  | 0.01                        | 97.4         | 1.0     |
| K           | 0.29                 | 0.281                        | 0.0001  | 0.005                       | 97.0         | 2.5     |
| Mg          | 0.28                 | 0.294                        | 0.0004  | 0.002                       | 94.7         | 0.8     |
| Na          | 0.052                | 0.052                        | 0.002   | 0.001                       | 99.6         | 2.2     |
| Si          | 5.61                 | 6.05                         | 0.004   | 0.06                        | 92.2         | 1.3     |
| Ti          | 0.18                 | 0.173                        | 0.00008 | 0.002                       | 96.4         | 1.4     |

Dilution factor: 513  
n = 10

# Precision and Accuracy in Coal Ash

19

| VS 7177-95 |                                   |                              |         |                             |              |         |
|------------|-----------------------------------|------------------------------|---------|-----------------------------|--------------|---------|
| Elements   | Certified values (%) <sup>*</sup> | Experimental values (%) n=10 | MDL (%) | Precision $t_{0,975;9}$ (%) | Accuracy (%) | RSD (%) |
| Al         | 14.33                             | 14.0                         | 0.01    | 0.1                         | 97.7         | 1.3     |
| Ca         | 3.49                              | 3.70                         | 0.02    | 0.06                        | 94.0         | 2.4     |
| Fe         | 3.83                              | 3.73                         | 0.003   | 0.03                        | 97.3         | 1.0     |
| K          | 0.49                              | 0.446                        | 0.0004  | 0.006                       | 91.0         | 1.9     |
| Mg         | 0.893                             | 0.874                        | 0.001   | 0.009                       | 97.9         | 1.5     |
| Na         | 0.10                              | 0.11                         | 0.01    | 0.01                        | 97.6         | 3.1     |
| Si         | 27.43                             | 28.6                         | 0.014   | 0.4                         | 95.7         | 2.1     |
| Ti         | 0.36                              | 0.357                        | 0.0003  | 0.003                       | 99.3         | 1.2     |

<sup>\*</sup>Measured as Al but reported as Al<sub>2</sub>O<sub>3</sub> on the certificate of analysis provided with material

Dilution factor: 1992  
n = 10

# Recovery for the 3 CRMs

20

| Elements | EOP 12-1-02<br>(brown coal fly ash) |               |          | NCS FC28127<br>(coal ash) |               |          | VS 7177-95<br>(coal)  |               |          |
|----------|-------------------------------------|---------------|----------|---------------------------|---------------|----------|-----------------------|---------------|----------|
|          | Spiked values<br>mg/L               | Recovery<br>% | RSD<br>% | Spiked values<br>mg/L     | Recovery<br>% | RSD<br>% | Spiked values<br>mg/L | Recovery<br>% | RSD<br>% |
| Ba       | 1                                   | 102           | 0.4      | 1                         | 102           | 1        | 1                     | 103           | 0.2      |
| K        | 3                                   | 91            | 3        | 5                         | 102           | 2        | 3                     | 91            | 4        |
| Mn       | 1                                   | 102           | 1        | 2                         | 100           | 0.7      | 1                     | 104           | 0.6      |
| P        | 1                                   | 99            | 1        | 1                         | 103           | 1        | 1                     | 101           | 3        |
| S        | 1                                   | 97            | 3        | 5                         | 100           | 1        | 1                     | 96            | 3        |
| Sr       | 1                                   | 96            | 1        | 1                         | 100           | 0.9      | 1                     | 98            | 3        |
| Zn       | 1                                   | 100           | 0.7      | 1                         | 101           | 0.9      | 1                     | 101           | 1        |

n = 5

## Conclusion

- Great overall accuracy:
  - Minimum (smaller value) 91.0% for K
  - Average  $96.6 \pm 2.3\%$
- Great overall recovery:
  - Minimum 91% for K
  - Average  $99.5 \pm 3.6\%$
- It is appropriate to use LeNeo instrument to prepare coal samples for ICP-OES analysis in accordance with industry norms



Questions?

[maprovencher@claisse.com](mailto:maprovencher@claisse.com)

Booth #722

claisse.com